InSpectra™ StO₂

SPOT CHECK

Model 300 Instructions for Use

ENGLISH

FRANÇAIS

DEUTSCH

ITALIANO

ESPAÑOL

Hutchinson
TECHNOLOGY™

Advancing Better Therapies

InSpectra™ StO₂

SPOT CHECK

Model 300 Instructions for Use

ENGLISH

FRANÇAIS

DEUTSCH

ITALIANO

ESPAÑOL

**Hutchinson
TECHNOLOGY™**

Advancing Better Therapies

InSpectra™ StO₂ Spot Check

Model Number: 300

Record:

Serial Number: _____

Purchase Date: _____

Hutchinson Technology Service Center

Hutchinson Technology Inc.
BioMeasurement Division
40 West Highland Park Drive NE
Hutchinson, MN 55350-9784
USA

Telephone: +1.800.419.1007 or +1.320.587.1926

Fax: +1.320.587.1555

Email: InSpectraStO2@hti.htch.com

www.htbiomeasurement.com

US Patents: 5879294, 6487343, 6473632, 6892006, 7460897, 7596397; others pending.

Copyright 2014. Hutchinson Technology Incorporated. Printed in the USA. All Rights Reserved.

InSpectra™ is a trademark of Hutchinson Technology Incorporated.

InSpectra™ is a registered trademark of Hutchinson Technology Inc. in the United States of America, the European Community, Canada, China, Japan, Saudi Arabia, Israel, New Zealand, Oman, Qatar, Australia, Lebanon and Russia and may be a registered trademark of Hutchinson Technology Inc. in other countries.

Table of Contents

InSpectra™ StO ₂ Spot Check	4
Introduction.....	4
Intended Use/Indications for Use	4
Contraindications	4
Safety.....	4
Safety Notifications.....	5
Rx Only.....	5
Symbols	6
How To Use.....	7
Component Set-Up.....	7
Features of the InSpectra™ StO ₂ Spot Check System.....	8
Battery	9
Charging the Battery.....	9
InSpectra™ StO ₂ Cable, model 3003	10
InSpectra™ StO ₂ Thenar Clip, model 1315	10
Power On/Off	11
Calibration.....	11
System Check.....	11
Patient Measurement.....	12
Positioning the InSpectra™ StO ₂ Clip	12
Data Display.....	13
Removing the InSpectra™ StO ₂ Clip	14
Maintenance	14
Cleaning the System.....	14
Software Updates.....	14
Maintenance.....	14
Disposal	15
Troubleshooting.....	15
Error Messages.....	15
Troubleshooting.....	16
Replacement Parts and Accessories	17
Technology.....	18
Technology Overview.....	18
Tissue Hemoglobin Index (THI).....	19
References.....	19
Technical Specifications	20
Warranty.....	22
Technical Addendum	23
Notes.....	28
Index	30

Introduction

The InSpectra™ StO₂ Spot Check uses near infrared spectroscopy to measure local oxygen saturation of hemoglobin in tissue (StO₂) and an index of total hemoglobin (THI). The absorption spectrum of light returned from a tissue sample varies primarily with oxyhemoglobin and deoxyhemoglobin concentration. The percent of StO₂ is a measure of hemoglobin oxygen saturation of blood contained in the volume of tissue illuminated by the near infrared light.

Tissue Hemoglobin Index (THI) is a quantified value that corresponds to the amount of hemoglobin present within the volume of tissue illuminated by the optical signal. Because the volume of tissue illuminated is comprised of extravascular and vascular tissue, THI measurements do not reflect blood hemoglobin content alone. The amount of tissue hemoglobin is influenced by blood hemoglobin concentration and microvasculature volume. With StO₂, THI serves as a hemoglobin signal strength indicator at the measurement site.

The InSpectra StO₂ Spot Check uses tissue absorbance values at different wavelengths. A reference measurement and scaling of the second derivative absorbance spectrum provide a means of internal calibration. Therefore, the InSpectra StO₂ Spot Check does not require recalibration before each use.

Intended Use/Indications for Use

Hutchinson Technology Incorporated's InSpectra™ StO₂ Spot Check is intended for use as a non-invasive monitoring system that measures an approximated value of percent hemoglobin oxygen saturation in the thenar skeletal muscle tissue (StO₂).

The InSpectra StO₂ Spot Check is indicated for use in monitoring patients during circulatory or perfusion examinations or when there is a suspicion of compromised circulation of the thenar skeletal muscle.

Contraindications

No contraindications are known.

Safety

Carefully read and review this entire manual to fully understand the device, its operation, and the potential risks to user and patient before use of the InSpectra™ StO₂ Spot Check System. Illustrations shown in this manual are intended to familiarize the user with the device. The StO₂ and THI values shown are intended solely as examples.

Safety Notifications

-  WARNINGS alert the user to potential serious outcomes (death, injury, or adverse events) to the patient or user and contain information to avoid injury.
-  CAUTIONS indicate conditions that may lead to malfunction or failure of the InSpectra™ StO₂ Spot Check System.
-  SAFETY NOTIFICATIONS help to:
- Identify the hazard.
 - Avoid the hazard.
 - Recognize the consequences.

Rx Only


Caution: Federal Law (USA) restricts this device to sale by, or on the order of, a physician.

Symbols

Safety Symbols

Rx Only	CAUTION: Federal Law (USA) restricts this device to sale by, or on the order of, a physician.		Follow instructions for use.
	Caution, consult accompanying documents.		Type BF Equipment – Defibrillator Proof
IPX0	Not protected against the ingress of fluid.		Class II Equipment
	Latex Free		Polarity
	Direct Current		WEEE Compliant

Device Symbols

	Reference Number		Serial Number
	Lot Number		
	Rechargeable Battery		Recycle
	Temperature Limitations	15% - 90% 60kPa - 106kPa	Relative Humidity Pressure
	Power Button		System Check

Agency Symbols

			Product Manufacturer
			Electrical Certification Conforms to UL STD 60601-1 Certified to CAN/CSA STD C22.2 NO. 601.1

Component Set-Up

Confirm all components below are present. If any component is missing or damaged, contact Hutchinson Technology Service Center.

Main components of the InSpectra™ StO₂ Spot Check System include:

- 1 InSpectra™ StO₂ Spot Check
- 2 InSpectra™ StO₂ Cable
- 3 InSpectra™ StO₂ Clip
- 4 Battery (rechargeable)
- 5 Charging Station
 - Power Supply
 - InSpectra™ StO₂ Spot Check model 300 Instructions for Use

The InSpectra™ StO₂ Spot Check can be used:

- as a hand-held portable device.
- in the Charging Station.
- mounted on an IV pole or flat vertical surface using the Mounting System accessory.

figure 7.1

WARNING:

- To ensure patient safety carefully route cabling to reduce the possibility of patient entanglement or strangulation.
- Do not place the InSpectra StO₂ Spot Check in any position that might allow it to fall on the patient.
- Explosion hazard. Do not use the InSpectra StO₂ Spot Check in the presence of flammable anesthetics or gases.

CAUTION:

- Do not use in the presence of ionizing radiations such as computed tomography (CT) or X-ray. Exposure to ionizing radiation will cause cumulative damage to the InSpectra StO₂ Spot Check and Cable resulting in inaccurate InSpectra StO₂ System Readings.
- The use of clips, cables, power supplies, and accessories other than those provided by Hutchinson Technology Inc. may result in damage and inaccurate readings.

Features of the InSpectra™ StO₂ Spot Check System

figure 8.1

figure 8.2

1 InSpectra™ StO ₂ Thenar Clip, model 1315	Reusable and easily cleaned Clip.
2 InSpectra StO ₂ Cable, model 3003	Optical cable features a push-button connection to the InSpectra StO ₂ Spot Check and an ejection button to easily release Clip.
3 Battery Indicator	Indicates the charge status of the battery. When a low battery is detected, the icon turns red and blinks. Indicator appears only when the InSpectra StO ₂ Spot Check is powered on.
4 InSpectra StO ₂ System Measurement	Two-digit display of percent tissue oxygen saturation (StO ₂) and THI; updated every two seconds.
5 THI Value (Tissue Hemoglobin Index)	
6 THI Signal Bars (Tissue Hemoglobin Index)	Provides a visual indicator of the THI signal strength.
7 Power Button	Used to Power On, Power Off and bring out of Sleep Mode.
8 System Check Button	When pressed, initiates the System Check to verify the function of the InSpectra StO ₂ Spot Check.
9 Battery	When new, the rechargeable lithium ion battery allows over 8 hours of continuous use.
10 Charging Station	Includes two sets of LED charge indicators to show status of two simultaneously charging batteries.

Battery

To install the battery in the InSpectra StO₂ Spot Check:

1. Place the battery in the back of the device as shown in Figure 9.1.
2. Slide battery in to engage.

NOTE: The battery provides a minimum of 8 hours of use when new a fully charged.

To remove:

1. Press button on the back of the battery.
2. Slide battery out to remove from device.

figure 9.1

Charging the Battery

To Set Up Charging Station:

1. Place on a flat surface next to the patient or mounted on an IV pole or flat vertical surface using the Mounting System accessory.

NOTE: Refer to the Mounting System Instructions for Use for additional information.

2. Using the power supply provided, plug into the receptacle on the side of the Charging Station and wall receptacle.

WARNING:

- Possible electrical shock hazard:
 - Do not connect to a wall receptacle if the power supply with the InSpectra StO₂ Spot Check is damaged.
 - Do not simultaneously touch a patient and charging contacts of the Charging Station.
- Possible Tripping Hazard. Do not allow the Power Supply to lie on the floor where it may become damaged or pose a tripping hazard.
- Do not position device so that it is difficult to disconnect from the AC outlet, as the disconnection of the power cord from the AC outlet is the method of separating the device from AC mains.

- 6 **CAUTION:** Use only the power supply provided with the InSpectra StO₂ Spot Check System to avoid damage to electrical components.

To Charge the Battery:

1. While installed in the InSpectra StO₂ Spot Check:
 - Place the InSpectra StO₂ Spot Check in the Charging Station as shown in Figure 9.2.

NOTE: Ensure the contacts protruding from the Charging Station are in contact with the device.

2. As an Independent Component:
 - Insert the battery into the charging station with the contacts pointing downward and facing out as shown in Figure 9.3.

When the battery properly inserted, the spare battery charge LEDs will indicate the battery status:

- Green – indicates the battery is fully charged.

figure 9.2

figure 9.3

Charging the Battery (continued)

- Yellow – indicates the battery is charging.
- Red – a battery charging error has occurred.
- Red and Yellow - the battery has overheated.

6 CAUTION: Possible battery discharge may result from extended non-use periods or storage. When returned to service, the InSpectra StO₂ Spot Check may require an extended charge time.

NOTE:

- InSpectra StO₂ Measurements can be taken while the InSpectra StO₂ Spot Check is in Charging Station.
- To ensure maximum battery runtime, fully charge the battery:
 - Upon receipt.
 - If the battery has not been fully charged in the past 30 days.
 - When the battery charge indicator turns red and begins to blink.
- It is recommended that the battery pack be stored and charged between 0°C–30°C with relative humidity less than 85%. If stored and charged at environmental conditions beyond these limits, overall battery capacity may be diminished, battery life may be shortened, and charge times may be lengthened.

InSpectra™ StO₂ Cable, model 3003

To Connect the InSpectra StO₂ Cable:

1. Align the button on the cable connector with the connector on the top of the InSpectra StO₂ Spot Check **1** as shown in Figure 10.1 and push together to connect.

To Remove InSpectra StO₂ Cable:

1. Press and hold the button on the cable connector while pulling the cable away from the InSpectra StO₂ Spot Check.

figure 10.1

6 CAUTION: Possible inaccurate measurements. Do not lift or carry the InSpectra™ StO₂ Spot Check by the Cable; doing so could damage the cable or device and result in inaccurate measurements.

InSpectra™ StO₂ Thenar Clip, model 1315

To Connect the InSpectra StO₂ Clip to the Cable:

1. Insert the two prongs on the InSpectra StO₂ Cable into the two holes on the top of the clip (**2**) as shown in Figure 10.2).
2. Press the cable and clip together to click into the locked position.

figure 10.2

To Remove InSpectra StO₂ Clip from Cable:

1. Pull back on the ejection lever.
2. Lift the Cable away from the Clip.

Power On/Off

To Power On InSpectra StO₂ Spot Check:

1. Press the power button on the front of the device.
2. The screen will illuminate and a series of white lines will gradually appear.
3. When complete, the battery indicator, StO₂ and THI labels will appear.

The InSpectra StO₂ Spot Check goes into Sleep Mode after 15 minutes to conserve battery power. Briefly press the power button to bring the InSpectra StO₂ Spot Check out of Sleep Mode.

To Power Off the InSpectra StO₂ Spot Check:

1. Press and hold the power button.
2. White lines will immediately appear on the screen and turn black as the device powers down, approximately 3 seconds.

➡ NOTE: Releasing the button at any time during shutdown will immediately restore power to the device.


figure 11.1

Calibration

The InSpectra™ StO₂ Spot Check System does not require user calibration. The InSpectra StO₂ Spot Check uses tissue absorbance values at different wavelengths. A reference measurement and scaling of the second derivative absorbance spectrum provide internal calibration.

System Check

The System Check provides a method to quickly check the ability of the InSpectra™ StO₂ Spot Check System to accurately calculate the InSpectra™ StO₂ System Measurement and THI.

Use the System Check:

- If the InSpectra StO₂ Spot Check System is exposed to ionizing radiation such as X-rays and computed tomography scanning.
- Upon replacement of the InSpectra™ StO₂ Cable or Clip.
- Whenever damage is suspected to the InSpectra StO₂ Spot Check or any of its components.

To use the System Check:

1. Position the InSpectra StO₂ Clip on the top of the InSpectra StO₂ Spot Check.
2. Press the System Check button.

➡ NOTE: If the InSpectra StO₂ Clip is not placed properly, an Error Screen will appear. Refer to the Error Message section for additional information.

figure 11.2

System Check (continued)

3. After 20 seconds, the InSpectra StO₂ System Measurement and THI will be displayed.

- The reading will appear in green if the System Check is within range.
- The reading will appear in red, if the InSpectra StO₂ Spot Check is out of range. Contact the Service Center.
- If the System Check reading failed, the fail icon will remain on screen. It will remain on screen until the reading passes.

6 CAUTION: If a failed System Check reading is obtained, do not measure a patient's StO₂. Refer to the Troubleshooting Guide.

4. All System Check readings will clear from the InSpectra StO₂ Spot Check after five 5 seconds.

To cancel a System Check reading at any time, press the System Check button to return to the main screen.


figure 12.1

Positioning the InSpectra™ StO₂ Clip

To Position the InSpectra StO₂ Clip:

1. Clean and dry the base of the thumb on the palmar side of the hand.

 **WARNING:** The InSpectra StO₂ Clip is for external use only and should only be used on intact skin.

2. Grasp clip by tab to open (Figure 12.2). Slide Clip over the patient's thumb and remove pressure. Slowly release pressure on tabs to secure placement on the patient.

figure 12.2

figure 12.3

figure 12.4

3. Ensure the InSpectra StO₂ Clip is in close contact with the palmar side of the hand (Figure 12.3 and 12.4).

4. Check THI signal strength on device screen.

 **NOTE:**

- A THI greater than 5.0 indicates sufficient hemoglobin to obtain an adequate signal in most circumstances.
- After positioning the clip, allow InSpectra StO₂ System Measurement approximately 20 seconds to settle prior to using or recording a reading.

Positioning the InSpectra™ StO₂ Clip (continued)

- Record the StO₂ reading before removing the clip from the measurement site. The InSpectra StO₂ Measurement is not retained by the device.

WARNING:

- Do not leave the InSpectra StO₂ Clip on the patient's thenar site for longer than 15 minutes. If continuous monitoring is desired, place the InSpectra™ StO₂ Tissue Oxygenation Monitor, model 650 and the InSpectra™ StO₂ Sensor, model 1615 for trended patient StO₂ measurements.
- Remove the InSpectra StO₂ Spot Check System from the patient prior to magnetic resonance imaging (MRI) scanning.

6 CAUTION:

- Do not use in the presence of ionizing radiations such as computed tomography (CT) or X-ray. Exposure to ionizing radiation will cause cumulative damage to the device resulting in inaccurately low InSpectra StO₂ System Readings.
- Ensure the InSpectra StO₂ Clip is fully attached to the patient to avoid high ambient light conditions that may result in the device not being able to take a measurement or provide inaccurate measurements.
- Excessive patient motion can cause large, rapid changes in InSpectra StO₂ Measurements.
- Inaccurate measurements can result from:
 - Contaminating the Clip on the patient interface area.
 - Placing the InSpectra StO₂ Clip over a hematoma or bony prominence.
 - Indocyanine green and methylene blue dyes in the blood.
 - Methemoglobin, carboxy-hemoglobin and other dyshemoglobins in the blood.
 - Placing the InSpectra StO₂ Clip over tissue obscured by ink such as tattoos, markers, pens.
- If a THI measurement is displayed when the InSpectra StO₂ Clip is not on a patient, refer to the Troubleshooting Section.

Data Display

The InSpectra™ StO₂ Spot Check screen consists of two areas:

Status Bar (top of screen)

The status bar at the top of the screen displays the status of the battery.

- The four segments of the battery icon provide information on the remaining battery capacity.
- At approximately 30 minutes remaining, a low battery icon will appear and blink to alert the user to recharge or replace the battery.

Data Display (main image area)

The data displayed provides the following live data updated every two (2) seconds:

- % StO₂ value
- Tissue Hemoglobin Index (THI) value
- THI signal bars
- Error Messages
- System Check results


figure 13.1

Removing the InSpectra™ StO₂ Clip

1. Gently apply pressure to the clip tabs until fully open.
2. Slide the Clip away from the patient's hand.

Cleaning the System

InSpectra™ StO₂ Spot Check, Cable, Clip and Charging Station:

1. Remove visible foreign material by wiping the surfaces with 70–90% Isopropyl Alcohol until visually clean (approximately 15 seconds).
2. Allow to air dry prior to returning it to operation.

6 CAUTION:

- Possible damage to electrical components. Do not immerse, spray, or drip fluid on any part of the InSpectra StO₂ Spot Check or Charging Station.
- Erroneous StO₂ readings may occur if gamma sterilization, autoclave, and Ethylene Oxide sterilization is used on any part of the InSpectra StO₂ Spot Check.
- Do not clean with acetone, ethyl alcohol or sodium hypochlorite (bleach) as this will damage the surface.
- Avoid use of cleaners that might scratch or dissolve plastics surfaces. Test the cleaner on a hidden surface for possible deterioration of exterior finishes.

Software Updates

Hutchinson Technology Incorporated may provide periodic software updates in order to increase the usability of the InSpectra™ StO₂ Spot Check.

Maintenance

The InSpectra StO₂ Cable, model 3003 and Clip, model 1315 can be replaced on site.

All other maintenance must be done by a Hutchinson Technology Incorporated Authorized Service Technician.

 **WARNING:** Possible Electrical Shock Hazard. Do not remove the cover of the InSpectra StO₂ Spot Check or Charging Station.

 **NOTE:** No modification to this equipment is allowed.

Disposal

The InSpectra™ StO₂ SpotCheck, Battery, Charging Station and Power Supply contain the following electronic components that must be disposed of in accordance with local regulations:

- Lithium ion cell in the battery.

WARNING: Chemicals from a broken LCD display panel or backlight are toxic when ingested; use caution when handling the device with a broken LCD screen.

Dispose the InSpectra™ StO₂ Cable and Clip per institutional guidelines.

Error Messages

Error message will remain on screen until the error is resolved.

For unresolved errors, contact the Hutchinson Technology Service Center.

	Graphic Display	Situation	Action
#001		System Check Not Seated Properly	1. Verify proper clip placement. 2. Rerun a System Check by pressing the System Check button.
#002		Defective Cable	1. Disconnect and reconnect cable. 2. If needed, connect a new cable.
#003		Over Temperature	1. Power off 2. Allow system to cool to ambient temperature. 3. Power on.
#004		Communication Error	1. Power Off, then Power On.
#005		Low LED Output	1. Power Off, then Power On.
#xxx		Software Error	1. Power Off, then Power On.

Troubleshooting

For additional information, contact your Hutchinson Technology representative/distributor or the Hutchinson Technology Service Center listed below.

Hutchinson Technology Service Center

Telephone: +1.800.419.1007 or +1.320.587.1926

Fax: +1.320.587.1555

Email: InSpectraStO2@hti.htch.com

Observation	Action
The InSpectra™ StO ₂ Spot Check does not power on	Battery may be discharged. 1. Charge the battery or place the device in the Charging Station.
No data on LCD screen.	1. Verify connection of the InSpectra StO ₂ Cable. 2. Verify connection of the InSpectra StO ₂ Clip. 3. Perform a System Check to verify the device is working. 4. Power Off, then Power On.
A pinwheel icon remain spinning on the screen. 	1. Verify connection of the InSpectra StO ₂ Cable. 2. Power Off, then Power On.
No InSpectra StO ₂ Measurement: doubledashes displayed(--).	Clip may not be placed properly over tissue with adequate hemoglobin content. 1. Adjust placement of Clip using THI as a signal strength indicator. Excessive residue under Clip. 1. Clean clip and patient skin where clip will be attached. Insufficient optical signal. The Cable or Clip may be damaged. 1. Perform a System Check to verify the device is working. 2. Replace the Cable.
Double dashes (- -) are displayed for InSpectra™ StO ₂ Measurement; values are displayed for THI.	Perform a System Check to verify the device is working.
A failed System Check was obtained.	Clip may not be placed properly on the System Check. 1. Verify placement of the InSpectra StO ₂ Clip on the System Check and rerun. InSpectra StO ₂ Cable may be damaged. 1. Replace InSpectra StO ₂ Cable and rerun System Check.
InSpectra StO ₂ Cable will not connect to InSpectra StO ₂ Spot Check.	Verify pins in cable connector are not bent. 1. If bent, carefully straighten pins. 2. Reconnect cable to the InSpectra StO ₂ Spot Check.
An Error Message remains on screen.	1. Refer to Error Messages.

Replacement Parts and Accessories

Please contact your Hutchinson Technology representative/distributor or the Hutchinson Technology Service Center to order replacement parts or accessories.

Hutchinson Technology Service Center

Telephone: +1.800.419.1007 or +1.320.587.1926

Fax: +1.320.587.1555

Email: InSpectraStO2@hti.htch.com

Part Number	Description
 5024730	InSpectra™ StO ₂ Spot Check System including: InSpectra™ StO ₂ Spot Check, Cable, Clip, Battery and Charging Station, Power Supply and Instructions for Use
Replacement Parts	
 5025585	InSpectra™ StO ₂ Spot Check Monitor
 5024735	InSpectra™ StO ₂ Clip
 5024732	InSpectra™ StO ₂ Cable
 5024734	Rechargeable Battery
 5024733	InSpectra™ StO ₂ Spot Check System model 300 Instructions for Use
 5026262	Charging Station with Power Supply
Accessories	
 5024745	Mounting System

Technology Overview

The InSpectra StO₂ Spot Check makes use of the characteristic light absorption properties of hemoglobin in the near infrared wavelength range. Near infrared spectroscopy (NIRS) relies on the relative transparency of living tissue to light in the 650–1000nm wavelength range. InSpectra StO₂ Spot Check measures tissue optical attenuation values at 680, 720, 760, and 800nm.

The light in the InSpectra StO₂ Cable contains the four wavelengths of light used for the InSpectra™ StO₂ System Measurement. The maximum depth of the tissue volume sampled is estimated to equal the distance between the clip's send and receive fibers. Cui, Kumar, and Chance (1991) confirmed that the mean measurement depth into the tissue is half of the clip spacing. The InSpectra™ StO₂ Clip model 1315: 15mm is designed to measure the proper depth of the tissue sampled in the thenar eminence.

There are two (2) light points on the face of the clip which send and receive the signal to and from the patient's tissue. The comparison of the receive signal from the patient and the receive feedback signal within the device is processed into a second derivative attenuation spectrum using a fixed wavelength gap point difference calculation. The resultant second derivative attenuation spectrum is sensitive to deoxyhemoglobin and oxyhemoglobin absorption. The absorption spectrum of light returned from a tissue sample varies mainly with oxyhemoglobin and deoxyhemoglobin concentration. Other chromophores have less effect.

Figure 18.1 illustrates inverted second derivative absorbance curves for hemoglobin at varying oxygen saturations.

Percent InSpectra StO₂ Measurements calculate an approximated value of percent hemoglobin oxygen saturation of blood contained in the volume of tissue illuminated by near infrared light. Currently, no measurement standard exists for measuring tissue hemoglobin oxygen saturation. A second derivative technique is used in order to quantify tissue hemoglobin oxygenation (StO₂).


figure 18.1 Second derivative absorbance curves for hemoglobin at varying oxygen saturation.

The second derivative spectral measurements at 720nm and 760nm are converted into a scaled second derivative spectrum value (2D720) using the 760nm measurement as the denominator within the calculated ratio. The second derivative algorithm method is insensitive to large changes in optical scattering. An empirically derived calibration curve, which relates the scaled second derivative calculated value to the displayed percent InSpectra StO₂ Measurement, is stored within the device.

Although a first derivative technique can remove baseline shifts from optical density spectra, it cannot remove sloping or tilting due to changes in wavelength dependent scattering, either over time or upon movement of the measurement clip to different tissue sites. Matcher and Cooper (1984) showed that the second derivative preprocessing of tissue attenuation measurements corrects for shifting and sloping baselines, as a result of measuring the attenuation of light in a tissue scattering environment. Research by Myers et al (2005) demonstrated that scaling the second derivative absorbance spectrum provides a measurement that is further insensitive to total hemoglobin and optical path length

Tissue Hemoglobin Index (THI)

Tissue Hemoglobin Index (THI), as measured in a range of 2–99 by the InSpectra™ StO₂ Spot Check, is a quantified value that corresponds to the amount of hemoglobin present within the volume of tissue illuminated by the optical signal of the InSpectra™ StO₂ Spot Check System. Because the volume of tissue illuminated is comprised of extravascular and vascular tissue, THI measurements do not reflect blood hemoglobin content alone. The amount of tissue hemoglobin present is influenced by blood hemoglobin concentration and microvasculature volume.

A THI reading of 10, measured by the InSpectra StO₂ Spot Check, indicates twice the signal strength as a reading of 5.0. Along with regional tissue oxygen saturation (StO₂), THI provides a method for determining the amount of hemoglobin present in the illuminated volume, and can therefore serve as a hemoglobin signal strength indicator at the measurement site. A THI greater than 5.0 indicates sufficient hemoglobin to obtain an adequate signal in most circumstances.

References

1. Cui W, Kumar C, Chance B. (1991). Experimental study of migration depth for the photons measured at sample surface. *Proc SPIE*, 1431, 180–191.
2. Mancini D, Bolinger L, Li H, Kendrick K, Chance B, Wilson JR. (1994). Validation of near-infrared spectroscopy in humans. *J Appl Physiol*, 77, 2740–2747.
3. Matcher SJ, Cooper CE. (1984). Absolute quantification of deoxyhaemoglobin concentration in tissue near infrared spectroscopy. *Phys Med Biol*, 39, 1295–1312.
4. Merrick MF, Pardue HL. (1986). Evaluation of absorption and first- and second-derivative spectra for simultaneous quantification of bilirubin and hemoglobin. *Clin Chem*, 32(4), 598–602.
5. Myers D, Anderson L, Seifert R, Ortner J, Cooper CE, Beilman G, Mowlem, JD. (2005). Noninvasive method for measuring local hemoglobin oxygen saturation in tissue using wide gap second derivative near infrared spectroscopy. *J BioMed Opt*, 10(3), 034017/1–18.

Technical Specifications

StO ₂ (Tissue Oxygen Saturation)	Range		0% to 99%
	Resolution		1%
	Display Update Rate		2 seconds
	Accuracy*	70% to 99%	1.63 StO ₂ units
		0% to 70%	2.94 StO ₂ units
THI (Tissue Hemoglobin Index)	Range		2 to 99 total range, 1 to 20 typical
	Resolution		0.1 THI units
	Display Update Rate		2 seconds
	Accuracy		Correlation to in vitro blood hemoglobin $R^2 > 0.95$
	StO ₂ Crosstalk		20% THI change with 100% full scale StO ₂ change
Alarms	Audiovisual		None
	Error Conditions		Visual flags for: • System Check Not Seated Properly • Defective Cable • Over Temperature • Communication Error • Low LED Output
Self-Test	Self-test occurs at Power On of InSpectra™ StO ₂ Spot Check.		
Environmental Parameters	Operation	Temperature	+10°C to +40°C (+50°F to +104°F)
		Relative Humidity	15% to 95%, non-condensing
		Pressure	106kPa to 70kPa
	Transport and Storage	Temperature	-18°C to +60°C (0°F to +140°F)
		Relative Humidity	15% to 90%, non-condensing
		Pressure	106kPa to 60kPa
Power	AC Voltage Input		100 to 240 VAC; 50 to 60 Hz Note: Use only the power supplied provided with the InSpectra StO ₂ Spot Check System
	Battery		• Type: Lithium Ion • Level of Charge Indicator • Charging status LED Indicators • Charges when connected to Charging Station regardless of power state
Physical Attributes	Size		• 18cm x 8cm x 3cm (7.25in x 3in x 1.25in)
	Weight		• 400g (15 oz)
	Cable		• Length: 1 m (3 ft) • Optical connection between Device and Clip
Display	Data Display		• Full Color Display • Live digital StO₂ and THI • THI bars

*Determined by in vitro measurements in a volume of blood as compared to a CO-oximeter (root mean square error).

Technical Specifications (continued)

Product Classification	Medical Device • EU: Class II device under the Medical Device Directive (93/42/EEC) as amended by 2007/47/EC per Annex IX, Section 3.2, Rule 10 • US: Class II device per section 870.2700 of 21 CFR • Canada: "Guidance for the Risk Based Classification System," Class II per Rule 10, subrule 1 • Electrical: Class 1, Type B F Defibrillation Proof, device is designed for continuous operation • Laser: Class 1, LASER Product	
Standards	Electrical and Constructional Safety	IEC 60601-1 Edition 3 v2005 UL-60601-1 First Edition v2003 AAMI ES60601-1 v2005 CSA-C22.2 No. 60101-1 Edition 3 v2008
	Electromagnetic Compatibility	IEC 60601-1-2 Edition 2.1 v2001 consolidated with Amendment 1:2004
	Enclosure	CEI IEC 60529 Edition 2 v1989; Amendment 1:1999, IPX0
	Laser Safety	IEC 60825-1: Edition 1.2 v2001 Corrigendum 1:2002; Class 1
	Shipping Stock	ASTM D 4169 – V2005: DC13, Assurance Level 2
Electromagnetic Emissions and Immunity	This equipment has been tested and found to comply with the limits for medical devices to IEC 60601-1-2:2001. This testing shows the device provides reasonable protection against harmful interference in a typical medical installation. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to other devices or is negatively impacted by other devices, the user is encouraged to try to correct the interference by one or more of the following measures: • Reorient or relocate the devices. • Increase the separation between the devices. • Connect the equipment to an outlet receptacle on a different circuit. • Contact the Service Center.	
InSpectra™ StO ₂ Clip, model 1315	Size	One size; fits left or right thenar eminence
	Attachment to Patient	Clip placement on thenar eminence Intact skin only Latex free
	Send to receive optical spacing	15mm (0.6 in)
	Biocompatibility Standards	Cytotoxicity: ISO 10993-5: v2009 Irritation and Sensitization: ISO 10993-10: v2002; Amendment 1:2006

Warranty

Hutchinson Technology Incorporated (HTI) warrants the InSpectra™ StO₂ Spot Check and its components according to the following terms and conditions.

InSpectra™ StO₂ Spot Check model 300

HTI warrants the InSpectra™ StO₂ Spot Check, model 300 to be free from defects in workmanship and material under normal use and service as specified in the Instructions For Use for a period of one year from the date of shipment by HTI or by or on behalf of an authorized distributor.

InSpectra™ StO₂ Clip model 1315

HTI warrants the InSpectra™ StO₂ Clip, model 1315 to meet HTI specifications at the time of purchase and to be free from original defects in material and workmanship under normal use and service as described in the Instructions For Use for a period of three months from the date of shipment by HTI or by or on behalf of an authorized distributor.

InSpectra™ StO₂ Cable model 3003

HTI warrants the InSpectra™ StO₂ Cable, model 3003 to meet HTI specifications at the time of purchase and to be free from original defects in material and workmanship under normal use and service as described in the Instructions For Use for a period of six months from the date of shipment by HTI or by or on behalf of an authorized distributor.

Rechargeable Battery for the InSpectra™ StO₂ Spot Check

HTI warrants the Rechargeable Battery to meet HTI specifications at the time of purchase and to be free from original defects in material and workmanship under normal use and service as described in the Instructions For Use for a period of thirty days from the date of shipment by HTI or by or on behalf of an authorized distributor.

Charging Station for the InSpectra™ StO₂ Spot Check

HTI warrants the Charging Station to meet HTI specifications at the time of purchase and to be free from original defects in material and workmanship under normal use and service as described in the Instructions For Use for a period of one year from the date of shipment by HTI or by or on behalf of an authorized distributor.

THIS WARRANTY IS EXCLUSIVE, AND HTI MAKES NO OTHER WARRANTY, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. Without limiting the foregoing, this is not a life-saving product or a diagnostic device. Any oral or written statement that is inconsistent with this Warranty shall be of no force or effect.

This Warranty does not cover abuse or physical damage. ANY UNAUTHORIZED WORK IMMEDIATELY VOIDS THIS WARRANTY. To obtain information or assistance, contact Customer Service.

This Warranty is not transferable and extends only to purchases directly from HTI or through an authorized representative.

Conditions

All claims under this Warranty must be submitted to HTI in writing within thirty days. Customer must obtain a Return Authorization Number from HTI by calling Customer Service. Customer must prepay shipping charges and is liable for loss or damage during shipment to HTI.

Upon receipt of the returned goods, HTI shall determine to its satisfaction that they have not been repaired or altered by anyone other than HTI; that their condition was not the result of abuse, neglect, or accident; that they have not been operated under conditions other than normal use; that all prescribed periodic maintenance and services have been performed; and that customer has met all other terms and conditions of this warranty.

Remedy

The liability of HTI under this Warranty is limited to one of the following, as HTI may elect in its sole discretion: replacing (with substantially similar, new or refurbished equipment) or repairing, issuing credit for, or refunding the purchase price of, products and components returned to HTI, adjusted to reflect their age and use.

All parts or equipment replaced will become the property of HTI.

THIS IS THE EXCLUSIVE REMEDY FOR ANY BREACH OF WARRANTY. The sole purpose of this remedy is to provide the Customer with the replacement (with substantially similar, new or refurbished equipment) or repair of nonconforming Product, or a credit or refund of a portion of the purchase price for such Product in the manner provided herein. This exclusive remedy shall not be deemed to have failed of its essential purpose so long as HTI is willing and able to repair or replace the nonconforming Product or credit or refund a portion of the purchase price therefore in the manner provided herein.

Limitation of Damages

IN NO EVENT SHALL HTI BE LIABLE TO CUSTOMER OR ANY THIRD PARTY FOR ANY INDIRECT, SPECIAL, CONSEQUENTIAL OR INCIDENTAL DAMAGES, INCLUDING WITHOUT LIMITATION, LOST PROFITS ARISING OUT OF THE USE OF THE PRODUCTS COVERED BY THIS WARRANTY OR CAUSED BY ANY DEFECT, FAILURE OR MALFUNCTION OF THOSE PRODUCTS WHETHER FOR BREACH OF ANY WARRANTY, FOR NEGLIGENCE, ON THE BASIS OF STRICT LIABILITY, OR OTHERWISE.

Technical Addendum

WARNING:

- The use of accessories, sensors and cables other than those specified may result in increased emission and/or decreased immunity of the InSpectra™ StO₂ Spot Check.
- The InSpectra StO₂ Spot Check System should not be used adjacent to or stacked with other equipment as close proximity with other electronics could result in reciprocal interference. If adjacent use or stacking is necessary, the InSpectra StO₂ Spot Check should be observed to verify normal operation in the configuration in which it will be used.

Guidance and Manufacturer's Declaration — Electromagnetic Emissions

The InSpectra™ StO₂ Spot Check is intended for use in the electromagnetic environment specified below. The customer or the user of the InSpectra StO₂ Spot Check should assure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment – Guidance
RF emissions CISPR 11	Group 1	The InSpectra StO ₂ Spot Check uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. The InSpectra StO ₂ Spot Check is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flickeremissions IEC 61000-3-3	Complies	

Technical Addendum

Guidance and Manufacturer's Declaration — Electromagnetic Immunity

The InSpectra™ StO₂ Spot Check is intended for use in the electromagnetic environment specified below. The customer or the user of the InSpectra StO₂ Spot Check should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Electrostatic discharge(ESD)	±6 kV contact	±6 kV contact	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
IEC 61000-4-2	±8 kV air	±8 kV air	
Electrical fast transient/burst	±2 kV for power supply lines	±2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
IEC 61000-4-4	±1 kV for input/output lines	±1 kV for input/output lines	
Surge	±1 kV differential mode	±1 kV differential mode	Mains power quality should be that of a typical commercial or hospital environment.
IEC 61000-4-5	±2 kV common mode	±2kV/commonmode	
Voltage dips, short interruptions and voltage variations on power supply input lines	<5 % U _T (>95 % dip in U _T) for 0.5 cycle	<5 % U _T (>95 % dip in U _T) for 0.5 cycle	Mains power quality should be that of a typical commercial or hospital environment. If the user of the InSpectra StO ₂ Spot Check requires continued operation during power mains interruptions, it is recommended that the InSpectra StO ₂ Spot Check be powered from an uninterruptible power supply or a battery.
IEC 61000-4-11	40 % U _T (60 % dip in U _T) for 5 cycles	40 % U _T (60 % dip in U _T) for 5 cycles	
	70 % U _T (30 % dip in U _T) for 25 cycles	70 % U _T (30 % dip in U _T) for 25 cycles	
	<5 % U _T (>95 % dip in U _T) for 5 seconds	<5 % U _T (>95 % dip in U _T) for 5 seconds	
Powerfrequency (50/60 Hz magnetic field)	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
IEC 61000-4-8			

NOTE: U_T is the a.c. mains voltage prior to application of the test level.

Technical Addendum

Guidance and Manufacturer's Declaration — Electromagnetic Immunity

(continued)

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150kHz to 80MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the InSpectra StO₂ Spot Check, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter, as portable and mobile RF communications equipment can effect medical electrical equipment. Recommended separation distance $d = 1.2\sqrt{P} \quad 150\text{kHz to } 80\text{MHz}$ Recommended separation distance $d = 1.2\sqrt{P} \quad 80\text{MHz to } 800\text{MHz}$ $d = 2.3\sqrt{P} \quad 800\text{MHz to } 2.5\text{GHz}$ Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80MHz to 2.5GHz	3 V/m	

➡ NOTE: At 80 MHz and 800 MHz, the higher frequency range applies.

➡ NOTE: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Technical Addendum

Recommended Separation Distances Between Portable and Mobile RF Communications Equipment and the InSpectra™ StO₂ Spot Check

The InSpectra StO₂ Spot Check is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the InSpectra StO₂ Spot Check can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the InSpectra StO₂ Spot Check as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150kHz to 80MHz $d = 1.2\sqrt{P}$	80MHz to 800MHz $d = 1.2\sqrt{P}$	800MHz to 2.5GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

➡ NOTE: At 80 MHz and 800 MHz, the higher frequency range applies.

➡ NOTE: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Technical Addendum

Guidance and Manufacturer's Declaration — Accessory Specifications

Part Number	Item
5024734	Rechargeable Battery
5026262	Charging Station with Power Supply
Use of accessories other than those specified above may result in increased emission and/or decreased immunity of the InSpectra StO ₂ Spot Check.	

Notes

[illegible]

Notes

[illegible]

Index

A

Accessories	7, 17
Algorithm	18
Ambient Light	13

B

Battery	3, 6, 7, 8, 9, 10, 11, 13, 15, 16, 17, 20, 22, 24
Battery Charge Indicator	10
Battery Discharge	10

C

Cable	7, 8, 10, 11, 14, 15, 16, 17, 18, 20, 22, 23, 24
Calibration	4, 11, 18
Charging Station	7, 8, 9, 10, 14, 15, 16, 17, 20, 21, 22, 27
Cleaning	14
Clip	7, 8, 10, 12, 13, 14, 15, 16, 17, 18, 20, 21, 22, 25
Communication Error	15, 20
Computed Tomography	7, 11, 13
Contraindications	4

D

Data Display	13, 20
Defective Cable	15, 20
Deoxyhemoglobin	4, 18
Display	8, 12, 13, 15, 16, 18, 20
Disposal	15
Double Dashes	16
Dyshemoglobins	13

E

Ejection Button	8
Electrical Certification	6
Electromagnetic Emissions	21, 23
Electromagnetic Immunity	21, 25
Environmental Parameters	20
Error Message	11, 13, 15, 16

H

Hematoma	13
Hemoglobin	4, 8, 13, 16, 18, 19, 20

I

Indications for Use	4
Indicators	8, 20
Instructions for Use	7, 9, 17, 22
Indocyanine Green	13
InSpectra StO ₂ Cable	7, 8, 10, 11, 14, 15, 16, 17, 18, 22
InSpectra StO ₂ Clip	7, 8, 10, 11, 12, 13, 14, 15, 16, 17, 18, 21, 22
Intended Use	4

L

LCD Screen	15, 16
LED Indicator	20
Lithium Ion Battery	8
Low LED Output	15, 20

Index

M

Magnetic resonance imaging	13
Methemoglobin	13
Methylene Blue.....	13
Model Number.....	2
Mounting System Accessory.....	7, 9

N

Near Infrared Spectroscopy.....	4, 18, 19
---------------------------------	-----------

O

Optical Signal	4, 16, 19
Over Temperature	15, 20
Oxygen Saturation	4, 8, 18, 19, 20
Oxyhemoglobin	4, 18

P

Patient Motion	13
Physical Attributes.....	20
Polarity	6
Power.....	6, 7, 8, 9, 11, 15, 16, 17, 20
Power Off	8, 11, 15, 16
Power On.....	8, 11, 15, 16, 20
Power Supply.....	7, 9, 15, 17, 27
Product Classification.....	21

R

Rechargeable Battery.....	6, 17, 22, 27
Replacement Parts	17

Rx Only.....	5
--------------	---

S

Safety.....	4, 5, 6, 7, 21
Serial Number	2, 6
Service Center	2, 7, 12, 15, 16, 17, 21
Sleep Mode	8, 11
Symbols	6
System Check	6, 8, 11, 12, 13, 15, 16, 20
System Check Button	8, 11, 12, 15
System Check Not Seated Properly	20

T

Technical Specifications	20, 21
Thenar Eminence.....	12, 18, 21
THI.....	4, 8, 11, 12, 13, 16, 19, 20
THI Signal Bars.....	8, 13
Tissue Hemoglobin Index.....	4, 8, 13, 19, 20
Tissue Oxygen Saturation.....	8, 19
Troubleshooting.....	12, 13, 16

V

Voltage.....	20
--------------	----

W

Warranty.....	22
WEEE	6

X

X-Ray	7, 11, 13
-------------	-----------

The InSpectra™ StO₂ Spot Check is intended for use as a noninvasive monitoring system that measures an approximated value of percent hemoglobin oxygen saturation in the tarsal skeletal muscle tissue (StO₂). Visit us online at www.htbiomeasurement.com for our full contact information and Instructions for Use.

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Références

1. Cui W, Kumar C, Chance B. (1991). Experimental study of migration depth for the photons measured at sample surface. Proc SPIE, 1431, 180–191.
2. Mancini D, Bolinger L, Li H, Kendrick K, Chance B, Wilson JR. (1994). Validation of near-infrared spectroscopy in humans. J Appl Physiol, 77, 2740–2747.
3. Matcher SJ, Cooper CE. (1984). Absolute quantification of deoxyhaemoglobin concentration in tissue near infrared spectroscopy. Phys Med Biol, 39, 1295–1312.
4. Merrick MF, Pardue HL. (1986). Evaluation of absorption and first-and-second-derivatives spectra for simultaneous quantification of bilirubin and hemoglobin. Clin Chem, 32(4), 598–602.
5. Myers D, Anderson L, Seifert R, Ortner J, Cooper CE, Beilman G, Mowlem, JD. (2005). Noninvasive method for measuring local hemoglobin oxygen saturation in tissue using wide gap second derivative near infrared spectroscopy. J BioMed Opt, 10(3), 034017/1–18.

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Literatur

1. CuiW, KumarC, ChanceB. (1991). Experimental study of migration depth for the photons measured at sample surface. *Proc SPIE*, 1431, 180–191.
2. ManciniD, BolingerL, LiH, KendrickK, ChanceB, WilsonJR. (1994). Validation of near-infrared spectroscopy in humans. *J Appl Physiol*, 77, 2740–2747.
3. MatcherSJ, CooperCE. (1984). Absolute quantification of deoxyhaemoglobin concentration in tissue near infrared spectroscopy. *Phys Med Biol*, 39, 1295–1312.
4. MerrickMF, PardueHL. (1986). Evaluation of absorption and first- and second-derivative spectra for simultaneous quantification of bilirubin and hemoglobin. *Clin Chem*, 32(4), 598–602.
5. MyersD, AndersonL, SeifertR, OrtnerJ, CooperCE, BeilmanG, Mowlem, JD. (2005). Noninvasive method for measuring local hemoglobin oxygen saturation in tissue using wide gap second derivative near infrared spectroscopy. *J BioMed Opt*, 10(3), 034017/1–18.

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